

Instructions for the Non-Looping Multi-Event Monitor

MONITOR SET-UP:

1. Install an AAA type battery.
2. Clear memory. See **RESETTING THE MONITOR** section.

WHEN SYMPTONS OCCUR

1. Be as still as possible while recording.
2. Hold down **ORANGE** key until sound is heard.
3. Place **RIGHT** and **LEFT** thumbs on the metal electrodes according to the markings and hold them there until sound stops.
4. At the end of the recording, the sound stops and the **GREEN** ECG indicator flashes for 30 seconds at a rate corresponding to the number of events stored. (up to 4)

Number of flashes equals number of events stored – One Flash = 1 stored event; Two flashes – 2 events; Three flashes = 3 events; Four flashes = 4 events

TRANSMITTING:

1. Call 1-800-747-4455 as soon as possible after each event. A technician will help you through the process. For best transmission, please avoid cordless and cell phones unless it is an emergency.
2. Verify that the slider switch on the side of the device is set to **TEL** position. (By rear label of device)
3. Place phone mouthpiece the opening in the center of the devices green circle.
4. Press the **GRAY** key until a sound is emitted.
5. Transmission takes approximately 30 – 120 seconds.
6. End of transmission: Sound stops, snooze alarm starts and **GREEN** indicator light flashes for 30 seconds.
7. Pick up the phone and listen to the technician for further instructions.

RESETTING THE MONITOR:

1. **RESET** monitor by holding both **ORANGE** and **GRAY** buttons down at the same time.
2. **HOLD** until you hear 5 beeps. Monitor is now reset and ready to record again.
3. Removing the battery will not erase the memory.

