

Instructions for the Non-Looping Multi-Event Monitor

MONITOR SET-UP:

1. Clear memory. See **RESETTING THE MONITOR** section.

WHEN SYMPTOMS OCCUR

1. At first sign of a symptom, place the monitor's electrode feet evenly against the chest on the upper part of sternum.
2. Press the red **EVENT** button.
3. The yellow **MEM** light will turn to a solid ON status to indicate recording and tones will sound to let you know that the monitor is recording properly.
4. You must remain still and try to hold the monitor as still as possible until the tones stop, approximately 30 seconds.

TRANSMITTING:

1. Call 1-800-747-4455 as soon as possible after each event. A technician will help you through the process. For best transmission, please avoid cordless and cell phones unless it is an emergency.
2. Lay the monitor on a flat surface and place the mouthpiece of the telephone over the speaker port in the center of the **IMD** logo.
3. Press and release the grey **SEND** button to begin transmission. A series of tone will sound, and when they stop, pick up the phone and speak with the technician again.



RESETTING THE MONITOR:

Erase An Event

1. After the monitor has completed the transmission, you should pick up the telephone and speak with the technician to verify the transmission was received successfully
2. If the transmission was successful, no further action is required. You can simply record another set of events as described in **WHEN SYMPTOMS OCCUR**

Manual Erase Procedure

1. Only upon instruction from the technician, erase the memory by holding down the **EVENT** button, and then press the **SEND** button once.
2. Release both buttons simultaneously.
3. The recorder will beep an approval tone and the yellow **MEM** light will flash to verify that the memory has been erased.

*When the memory of the monitor is full, the yellow **MEM** light will flash 10 times at the end of the last event recording, and beep 3 times. If an attempt is made to record further events, an error message will sound accompanied by 3 flashes of the yellow **MEM** light*